

Inskam Find B WiFi Digital Otoscope for Ear Cleaning



Precautions

1. Do not use this tool if you have an active ear infection, ruptured eardrum or redness. Discontinue use immediately and seek medical care whenever.
2. Please read the instructions carefully. We do not assume responsibility for any injuries that might be caused by our products. If you have ear drainage, discharge, ear pain, ear perforation or injury, please consult with doctors before using.
3. The downloaded app helps you visualize and remove ear wax from the ear canal. Please keep the ear spoon at a safe distance from the eardrum.
4. Please be careful when using it to avoid getting too deep into the ear canal and causing injury to the ear canal.
5. If the ear wax is hard and stubborn to remove, please use ear drops for a few days to soften it before removing. Do not press too hard on your ear canal when removing the wax to avoid any scraping.
6. If there is a significant buildup of earwax when first using the device, it's best to remove the wax in stages to avoid irritation.
7. After the initial wax cleaning, use the device as required to avoid any further build-up of wax.
8. Always disinfect the accessories and device with alcohol after use. When cleaning the lens, please wipe carefully with a professional alcohol wipe cotton swab. It is forbidden to soak in water or wash the whole body.
9. It is forbidden to use this device while charging.
10. It is prohibited for children to use this product alone to avoid careless harm.

Steps to use on the phone

Step 1

For IOS:

Download and install the App “Soulear” from the App Store.

For Android:

Download and install the App “Soulear” from Google Play Store.



Optional:

Scan the QR Code below to download the App.


Step 2

Long press the switch key to open the device, after 5 seconds, then step into the next operation.




Step 3

Open the "Soulear" App, click the "Connect Device", follow the pop-up prompt and click to connect now to jump to the mobile WIFI setting interface.




Step 4

Find "Soulear-xxxx" on the WIFI setting interface and click to connect.



Step 5


After the mobile phone displays "Soulear-xxxx" WIFI is connected, return to the App again and click the "Connected equipment" to enter the device.



After connecting to the product WIFI connection, you will be asked some questions:

IOS:


"Soulear" would like to find and connect to devices on your local network. Suggestion: Click "OK"



Allow "Soulear "to use wireless data?

Suggestion:

Click "WLAN&Cellular/WLAN Only"




Android:

If your mobile phone displays the following prompts, WiFi has no data connection, whether to continue trying to use the WiFi? Whether to keep the connection? Whether to switch WiFi?

Suggestion:

Click the “Keep”/”CONNECT”/”Decline” option.



Troubleshooting

Cannot store photos or videos on your phone?

- Let Soulear permit the authority of photos, media and files in your device.

Camera has no response, camera doesn't light up?

- Please charge the device.

App crashed?

- Please uninstall the App and scan the QR code to install.

Live image does not show after launching the app?

- Make sure you connect the correct WIFI ID "Soulear-xxxx".
- The question: This WLAN network has no Internet access. Connect anyway? Please be sure to choose "Connect". If you choose "CANCEL" at the first time, you need to forget the network of the "Soulear-xxxx" WIFI on your phone, then connect with this WIFI again. Continue to open the app and choose "CONNECT" of this question, so that it can work properly.
- Enter your mobile device setting and allow the "Soulear" App to access "WLAN" or "WLAN & Cellular Data".
- Close the App and open it again.
- Turn the mobile phone data on, then off again.

Ear endoscope's red light on?

- Please charge the product.

WIFI connection break?

- Please restart the ear endoscope, open the App, connect to WIFI.

When iPhone prompts "cannot join this WIFI, the network is temporarily unavailable, please try again later" ?

- When the above situation occurs, please try to turn off the product, turn it on again after 30s, and reconnect to WIFI to connect and use.

When an Android phone connects to WIFI, it prompts "Network access is full" ?

- When the above situation occurs, please try to turn off the product, turn it on again after 30s, and reconnect to WIFI to connect and use.

The power indicator of the ear pick shows red?

- Please use the matching charging cable to charge.